

## 兔抗 MYH14 多克隆抗体

中文名称: 兔抗 MYH14 多克隆抗体

英文名称: Anti-MYH14 rabbit polyclonal antibody

别名: MRIP; RIP3; M-RIP; RHOIP3; p116Rip

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: MYH14

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and capping. By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MYH14
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MYH14.
<b>Full name:</b>	myosin phosphatase Rho interacting protein
<b>Synonyms :</b>	MRIP; RIP3; M-RIP; RHOIP3; p116Rip
<b>SwissProt:</b>	Q7Z406
<b>WB Predicted band size:</b>	228 kDa

<b>WB Positive control:</b>	293 cells lysate
<b>WB Recommended dilution:</b>	500-3000

