

兔抗 MYH4 多克隆抗体

中文名称: 兔抗 MYH4 多克隆抗体

英文名称: Anti-MYH4 rabbit polyclonal antibody

别名: MYH2B; MyHC-2B; MyHC-IIb

储存: 冷冻 (-20℃) 避光

抗原: MYH4

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Muscle contraction.
Applications:	IHC
Name of antibody:	MYH4
Immunogen:	Synthesized peptide derived from C-terminal of human MYH4.
Full name:	myosin, heavy chain 4, skeletal muscle
Synonyms :	MYH2B; MyHC-2B; MyHC-IIb
SwissProt:	Q9Y623
IHC positive control:	Human skeletal muscle tissue
IHC Recommend dilution:	50-100

