

兔抗 MYL3 多克隆抗体

中文名称：兔抗 MYL3 多克隆抗体

英文名称： Anti-MYL3 rabbit polyclonal antibody

别名： CMH8; VLC1; MLC1V; MLC1SB

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MYL3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	MYL3 encodes myosin light chain 3, an alkali light chain also referred to in the literature as both the ventricular isoform and the slow skeletal muscle isoform. Mutations in MYL3 have been identified as a cause of mid-left ventricular chamber type hypertrophic cardiomyopathy.
Applications:	ELISA, WB, IHC
Name of antibody:	MYL3
Immunogen:	Fusion protein of human MYL3
Full name:	myosin, light chain 3, alkali; ventricular, skeletal, slow
Synonyms :	CMH8; VLC1; MLC1V; MLC1SB
SwissProt:	P08590

ELISA Recommended dilution:	2000-10000
IHC positive control:	Human lung cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	22 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	1000-5000



