

兔抗 MYL6 多克隆抗体

中文名称：兔抗 MYL6 多克隆抗体

英文名称： Anti-MYL6 rabbit polyclonal antibody

别名： LC17; ESMLC; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM; LC17-GI; LC17-NM

抗原： MYL6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain that is expressed in smooth muscle and non-muscle tissues. Genomic sequences representing several pseudogenes have been described and two transcript variants encoding different isoforms have been identified for this gene.
Applications:	ELISA, WB, IHC

Name of antibody:	MYL6
Immunogen:	Fusion protein of human MYL6
Full name:	myosin, light chain 6, alkali, smooth muscle and non-mu scl
Synonyms :	LC17; ESMLC; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM; LC17-GI; LC17-NM
SwissProt:	P60660
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	17 kDa
WB Positive control:	Mouse stomach tissue
WB Recommended dilution:	200-1000



