

兔抗 MYL9 多克隆抗体

中文名称：兔抗 MYL9 多克隆抗体

英文名称： Anti-MYL9 rabbit polyclonal antibody

别名： LC20, MLC2, MRLC1, MYRL2, MLC-2C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MYL9

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Myosin, a structural component of muscle, consists of two heavy chains and four light chains. The protein encoded by this gene is a myosin light chain that may regulate muscle contraction by modulating the ATPase activity of myosin heads. The encoded protein binds calcium and is activated by myosin light chain kinase. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	MYL9
Immunogen:	Synthetic peptide of human MYL9
Full name:	myosin, light chain 9, regulatory

Synonyms :	LC20, MLC2, MRLC1, MYRL2, MLC-2C
SwissProt:	P24844
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human liver cancer and Human cervical cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	14 kDa
WB Positive control:	Jurkat and 231 cells
WB Recommended dilution:	1000-5000



