

兔抗 MYO1D 多克隆抗体

中文名称: 兔抗 MYO1D 多克隆抗体

英文名称: Anti-MYO1D rabbit polyclonal antibody

别名: myr4; PPP1R108

储存: 冷冻 (-20°C) 避光

抗原: MYO1D

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Myosins are actin-based motor molecules with ATPase activity . Unconventional myosins serve in intracellular movements. Their highly divergent tails are presumed to bind to membranes and organelles, which would be moved relative to actin filaments. By similarity.
Applications:	IHC
Name of antibody:	MYO1D
Immunogen:	Synthesized peptide derived from internal of human MYO1D.
Full name:	myosin ID
Synonyms :	myr4; PPP1R108
SwissProt:	O94832

IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

