

兔抗 MYOD1 (Phospho-Ser200)多克隆抗体

中文名称：兔抗 MYOD1 (Phospho-Ser200)多克隆抗体

英文名称：Anti-MYOD1 (Phospho-Ser200) rabbit polyclonal antibody

别名：PUM; MYF3; MYOD; bHLHc1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：MYOD1 (Phospho-Ser200)

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	MyoD encodes a nuclear protein that belongs to the basic helix-loop-helix family of transcription factors and the myogenic factors subfamily. It regulates muscle cell differentiation by inducing cell cycle arrest, a prerequisite for myogenic initiation. The protein is also involved in muscle regeneration. It activates its own transcription which may stabilize commitment to myogenesis.
Applications:	WB
Name of antibody:	MYOD1 (Phospho-Ser200)
Immunogen:	Synthetic peptide of human MYOD1 (Phospho-Ser200)
Full name:	myogenic differentiation 1 (Phospho-Ser200)

Synonyms :	PUM; MYF3; MYOD; bHLHc1
SwissProt:	P15172
WB Predicted band size:	40 kDa
WB Positive control:	Hela cells untreated or treated with EGF
WB Recommended dilution:	500-1000

