

## 绵羊抗 FITC 多克隆抗体

中文名称： 绵羊抗 FITC 多克隆抗体

英文名称： Anti-FITC sheep polyclonal antibody

相关类别： 标签抗体

储 存： 冷冻（-20℃）

宿 主： Sheep

抗 原： FITC

反应种属： FITC

标 记 物： Unconjugate

克隆类型： sheep polyclonal

### 技术规格

<b>Background:</b>	Fluorescein isothiocyanate (FITC) is a derivative of fluorescein used in wide-ranging applications including flow cytometry. FITC is the original fluorescein molecule functionalized with an isothiocyanate reactive group (-N=C=S), replacing a hydrogen atom on the bottom ring of the structure. This derivative is reactive towards nucleophiles including amine and sulfhydryl groups on proteins.
<b>Applications:</b>	ELISA
<b>Name of antibody:</b>	FITC
<b>Immunogen:</b>	KLH conjugated FITC
<b>Full name:</b>	Fluorescein isothiocyanate
<b>CAS No.:</b>	27072-45-3



<b>ELISA Recommended dilution:</b>	2000-4000
------------------------------------	-----------