

兔抗 GIGYF1 多克隆抗体

- 中文名称：兔抗 GIGYF1 多克隆抗体
- 英文名称：Anti-GIGYF1 rabbit polyclonal antibody
- 别名：GRB10 interacting GYF protein 1; GYF1; PERQ1
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 抗原：GIGYF1
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the gyf family of adaptor proteins. The encoded protein contains a gyf protein interaction domain. It binds growth factor receptor bound 10, another adaptor protein that binds activated insulin-like growth factor 1 and insulin receptors and regulates receptor signaling.
Applications:	ELISA, IHC
Name of antibody:	GIGYF1
Immunogen:	Synthetic peptide of human GIGYF1
Full name:	GRB10 interacting GYF protein 1
Synonyms:	GYF1; PERQ1

SwissProt:	O75420
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	50-100

