

兔抗 PXN(Phospho-Ser178) 多克隆抗体

中文名称：兔抗 PXN(Phospho-Ser178) 多克隆抗体

英文名称： Anti-PXN(Phospho-Ser178) rabbit polyclonal antibody

别名： PXN

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PXN(Phospho-Ser178)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011]
Applications:	WB

Name of antibody:	PXN(Phospho-Ser178)
Immunogen:	Peptide sequence around phosphorylation site of Serine 178(A-L-S(p)-P-L) derived from Human Paxillin.
Full name:	paxillin
Synonyms :	PXN
SwissProt:	P49023
WB Predicted band size:	65 kDa
WB Positive control:	HT29 cells lysates
WB Recommended dilution:	500-1000

