

FITC 标记链霉亲和素

中文名称： FITC 标记链霉亲和素

英文名称： FITC-conjugated Streptavidin

相关类别： 抗原

储 存： 冷冻（-20℃）

概述

FITC-conjugated Streptavidin

技术规格

Full name	Streptavidin
Applications	1/100-200 for FCM
Background	Streptavidin is a tetrameric protein composed of identical subunits. The preparation contains an N- and C-terminal shortened variant (core streptavidin) with improved properties concerning homogeneity, solubility, resistance towards proteolytic degradation and accessibility of the biotin binding pocket as compared to native streptavidin. The high affinity recognition of biotin and biotinylated molecules has made streptavidin one of the most important components in diagnostics and laboratory kits.