

## 兔抗 GIMAP4 多克隆抗体

中文名称：兔抗 GIMAP4 多克隆抗体

英文名称： Anti-GIMAP4 rabbit polyclonal antibody

别名： IAN1; IAN-1; IMAP4; MSTP062

抗原： GIMAP4

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. The encoded protein of this gene may be negatively regulated by T-cell acute lymphocytic leukemia 1 (TAL1). In humans, the IAN subfamily genes are located in a cluster at 7q36.1.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	GIMAP4
<b>Immunogen:</b>	Synthesized peptide derived from internal of human GIMAP4.
<b>Full name:</b>	GTPase, IMAP family member 4

<b>Synonyms :</b>	IAN1; IAN-1; IMAP4; MSTP062
<b>SwissProt:</b>	Q9NUV9
<b>WB Predicted band size:</b>	38 kDa
<b>WB Positive control:</b>	Hela cells, HT-29 cells and RAW264.7 cells lysates
<b>WB Recommended dilution:</b>	500-3000

