

## 兔抗 GINS1 多克隆抗体

中文名称：兔抗 GINS1 多克隆抗体

英文名称：Anti-GINS1 rabbit polyclonal antibody

别名：PSF1; RP4-691N24.2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GINS1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	DNA replication complex GINS protein PSF1 is a protein that in humans is encoded by the GINS1 gene. The GINS complex plays an essential role in the initiation of DNA replication, and progression of DNA replication forks. GINS complex seems to bind preferentially to single-stranded DNA. GINS1 is essential for function. Component of the GINS complex which is a heterotetramer of GINS1, GINS2, GINS3 and GINS4. Forms a stable subcomplex with GINS4. GINS complex interacts with DNA primase in vitro.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GINS1

<b>Immunogen:</b>	Synthetic peptide of human GINS1
<b>Full name:</b>	GINS complex subunit 1 (Psf1 homolog)
<b>Synonyms :</b>	PSF1; RP4-691N24.2
<b>SwissProt:</b>	Q14691
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

