

兔抗 GINS3 多克隆抗体

中文名称：兔抗 GINS3 多克隆抗体

英文名称： Anti-GINS3 rabbit polyclonal antibody

别名： GINS complex subunit 3; PSF3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GINS3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein subunit of the GINS heterotrimeric complex, which is essential for the initiation of DNA replication and replisome progression in eukaryotes. Alternatively spliced transcript variants encoding distinct isoforms have been described.
Applications:	ELISA, IHC
Name of antibody:	GINS3
Immunogen:	Fusion protein of human GINS3
Full name:	GINS complex subunit 3
Synonyms:	PSF3
SwissProt:	Q9BRX5

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	50-100



