

兔抗 N4BP2L2 多克隆抗体

中文名称: 兔抗 N4BP2L2 多克隆抗体

英文名称: Anti-N4BP2L2 rabbit polyclonal antibody

别名: CG005; PFAAP5; 92M18.3

抗原: N4BP2L2

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	3 isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	N4BP2L2
Immunogen:	Synthesized peptide derived from internal of human N4BP2L2.
Full name:	NEDD4 binding protein 2-like 2
Synonyms :	CG005; PFAAP5; 92M18.3
SwissProt:	Q92802
WB Predicted band size:	67 kDa
WB Positive control:	LOVO cells lysate

WB Recommended dilution: 500-3000

