

兔抗 PXN 多克隆抗体

中文名称：兔抗 PXN 多克隆抗体

英文名称： Anti-PXN rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： PXN

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445).
Applications:	ELISA, WB, IHC
Name of antibody:	PXN
Immunogen:	Synthetic peptide of human PXN
Full name:	paxillin
SwissProt:	P49023

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human tonsil
IHC Recommend dilution:	30-150
WB Predicted band size:	65 kDa
WB Positive control:	HT-29 cell lysate
WB Recommended dilution:	200-1000



