

兔抗 PYCR2 多克隆抗体

中文名称：兔抗 PYCR2 多克隆抗体

英文名称： Anti-PYCR2 rabbit polyclonal antibody

别名： pyrroline-5-carboxylate reductase 2; HLD10; P5CR2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PYCR2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the pyrroline-5-carboxylate reductase family. The encoded mitochondrial protein catalyzes the conversion of pyrroline-5-carboxylate to proline, which is the last step in proline biosynthesis. Alternatively spliced transcript variants have been described for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	PYCR2
Immunogen:	Fusion protein of human PYCR2

Full name:	pyrroline-5-carboxylate reductase 2
Synonyms:	HLD10; P5CR2
SwissProt:	Q96C36
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and Human tonsil
IHC Recommend dilution:	50-100
WB Predicted band size:	34 kDa
WB Positive control:	293T, HepG2 and K562 cell lysates
WB Recommended dilution:	500-2000



