

兔抗 GIT1 多克隆抗体

中文名称：兔抗 GIT1 多克隆抗体

英文名称： Anti-GIT1 rabbit polyclonal antibody

抗 原： GIT1

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	GTPase-activating protein for the ADP ribosylation factor family. May serve as a scaffold to bring together molecules to form signaling modules controlling vesicle trafficking, adhesion and cytoskeletal organization. Increases the speed of cell migration, as well as the size and rate of formation of protrusions, possibly by targeting PAK1 to adhesions and the leading edge of lamellipodia. Sequesters inactive non-tyrosine-phosphorylated paxillin in cytoplasmic complexes.
Applications:	WB, IHC
Name of antibody:	GIT1
Immunogen:	Synthesized peptide derived from C-terminal of human GIT1.

Full name:	G protein-coupled receptor kinase interacting ArfGAP 1
Synonyms :	0
SwissProt:	Q9Y2X7
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	84 kDa
WB Positive control:	Hela cells, Jurkat cells and HepG2 cells lysates
WB Recommended dilution:	500-3000



