

兔抗 PYGB 多克隆抗体

中文名称：兔抗 PYGB 多克隆抗体

英文名称：Anti-PYGB rabbit polyclonal antibody

别名：GPBB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PYGB

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a glycogen phosphorylase found predominantly in the brain. The encoded protein forms homodimers which can associate into homotetramers, the enzymatically active form of glycogen phosphorylase. The activity of this enzyme is positively regulated by AMP and negatively regulated by ATP, ADP, and glucose-6-phosphate. This enzyme catalyzes the rate-determining step in glycogen degradation
Applications:	ELISA, WB

Name of antibody:	PYGB
Immunogen:	Synthetic peptide of human PYGB
Full name:	Phosphorylase, glycogen; brain
Synonyms :	GPBB
SwissProt:	P11216
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	97 kDa
WB Positive control:	SKOV3 cells
WB Recommended dilution:	200-1000

