

## 兔抗 NABP1 多克隆抗体

中文名称：兔抗 NABP1 多克隆抗体

英文名称： Anti-NABP1 rabbit polyclonal antibody

别名： nucleic acid binding protein 1; SSB2; OBFC2A; SOSS-B2

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： NABP1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Single-stranded DNA (ssDNA)-binding proteins, such as OBFC2A, are ubiquitous and essential for a variety of DNA metabolic processes, including replication, recombination, and detection and repair of damage (Richard et al., 2008 [PubMed 18449195]).
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NABP1
<b>Immunogen:</b>	Fusion protein of human NABP1
<b>Full name:</b>	nucleic acid binding protein 1
<b>Synonyms:</b>	SSB2; OBFC2A; SOSS-B2
<b>SwissProt:</b>	Q96AH0

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	Jurkat cell, Mouse liver tissue lysates
<b>WB Recommended dilution:</b>	500-2000



