

## 兔抗 GJA3 多克隆抗体

中文名称：兔抗 GJA3 多克隆抗体

英文名称： Anti-GJA3 rabbit polyclonal antibody

别名： CX46; CZP3; CTRCT14

抗原： GJA3

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a connexin and is a component of lens fiber gap junctions. Defects in this gene are a cause of zonular pulverulent cataract type 3 (CZP3).
<b>Applications:</b>	IHC
<b>Name of antibody:</b>	GJA3
<b>Immunogen:</b>	Synthesized peptide derived from internal of human GJA3.
<b>Full name:</b>	gap junction protein, alpha 3, 46kDa
<b>Synonyms :</b>	CX46; CZP3; CTRCT14
<b>SwissProt:</b>	Q9Y6H8
<b>IHC positive control:</b>	Human lung carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100

