

兔抗 6X His 多克隆抗体

中文名称: 兔抗 6X His 多克隆抗体

英文名称: Anti-6×His rabbit polyclonal antibody

别 名: 6 His epitope tag

相关类别: 标签抗体

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: 6×His

反应种属: Fusion proteins containing the HIS Tag

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

	The 6 X histidine motif is used as a tag on many rec
	ombinant proteins to facilitate purification.It can then
	be purified using immobilized metal ions,traditionally
	nickel,in a procedure termed immobilized-metal affinit
Background:	y chromatography. Further analysis of the protein can
	be achieved using an antibody specific to the 6 X hist
	idine in techniques such as western blotting. The antib
	ody recognizes the His-tag fused to targeted proteins
	in transfected or transformed cells.
Applications:	WB
Name of antibody:	6×His



Immunogen:	Synthetic peptide: 18×His conjugated to KLH
Full name:	НННННН Tag
Synonyms:	6 His epitope tag
WB Predicted band size:	59 kDa
WB Positive control:	Fusion proteins containing the HIS Tag
WB Recommended dilution:	1000-5000

