

兔抗 NAIF1 多克隆抗体

- 中文名称: 兔抗 NAIF1 多克隆抗体
- 英文名称: Anti-NAIF1 rabbit polyclonal antibody
- 别名: nuclear apoptosis inducing factor 1; C9orf90; bA379C10.2
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 抗原: NAIF1
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	Induces apoptosis.
Applications:	ELISA, IHC
Name of antibody:	NAIF1
Immunogen:	Fusion protein of human NAIF1
Full name:	nuclear apoptosis inducing factor 1
Synonyms:	C9orf90; bA379C10.2
SwissProt:	Q69YI7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and Human prostate cancer

IHC Recommend dilution: 50-100



