

兔抗 RAB3C 多克隆抗体

中文名称：兔抗 RAB3C 多克隆抗体

英文名称：Anti-RAB3C rabbit polyclonal antibody

别名：RAB3C, member RAS oncogene family

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RAB3C

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the RAS oncogene family and encodes a small GTPase. Other similar small GTPases are known to be involved in vesicle trafficking, and the encoded protein was shown to play a role in recycling phagocytosed MHC class 1 complexes to the cell surface. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	RAB3C

Immunogen:	Synthetic peptide of human RAB3C
Full name:	RAB3C, member RAS oncogene family
SwissProt:	Q96E17
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	26 kDa
WB Positive control:	Human fetal brain tissue and 293T cells
WB Recommended dilution:	500-2000



