

兔抗 GJB4 多克隆抗体

中文名称: 兔抗 GJB4 多克隆抗体

英文名称: Anti-GJB4 rabbit polyclonal antibody

别名: EKV; CX30.3

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GJB4

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a transmembrane connexin protein that is a component of gap junctions. Mutations in this gene have been associated with erythrokeratoderma variabilis, progressive symmetric erythrokeratoderma and hearing impairment.
Applications:	ELISA, WB, IHC
Name of antibody:	GJB4
Immunogen:	Synthetic peptide of human GJB4
Full name:	gap junction protein, beta 4, 30.3kDa
Synonyms :	EKV; CX30.3
SwissProt:	Q9NTQ9

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human colon cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	30 kDa
WB Positive control:	Mouse skin tissue
WB Recommended dilution:	200-1000



