

兔抗 RAB3IP 多克隆抗体

中文名称：兔抗 RAB3IP 多克隆抗体

英文名称： Anti-RAB3IP rabbit polyclonal antibody

别名： RAB3A interacting protein; RABIN3; RABIN8

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RAB3IP

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Guanine nucleotide exchange factor (GEF) which may activate RAB8A and RAB8B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form. Mediates the release of GDP from RAB8A and RAB8B but not from RAB3A or RAB5. Modulates actin organization and promotes polarized transport of RAB8A-specific vesicles to the cell surface. Together with RAB11A, RAB8A, the exocyst complex, PARD3, PRKCI, ANXA2, CDC42 and DNMBP promotes

	transcytosis of PODXL to the apical membrane initiation sites (AMIS), apical surface formation and lumenogenesis.
Applications:	ELISA, IHC
Name of antibody:	RAB3IP
Immunogen:	Fusion protein of human RAB3IP
Full name:	RAB3A interacting protein
Synonyms:	RABIN3; RABIN8
SwissProt:	Q96QF0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-100

