

兔抗 RAB4A 多克隆抗体

中文名称：兔抗 RAB4A 多克隆抗体

英文名称：Anti-RAB4A rabbit polyclonal antibody

别名：RAB4; HRES1; HRES-1; HRES-1/RAB4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RAB4A

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the largest group in the Ras superfamily of small GTPases, which regulate membrane trafficking. The encoded protein is associated with early endosomes and is involved in their sorting and recycling. The protein also plays a role in regulating the recycling of receptors from endosomes to the plasma membrane. Alternatively spliced transcript variants have been observed for this gene.
Applications:	ELISA, IHC

Name of antibody:	RAB4A
Immunogen:	Full length fusion protein
Full name:	RAB4A, member RAS oncogene family
Synonyms:	RAB4; HRES1; HRES-1; HRES-1/RAB4
SwissProt:	P20338
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	30-150

