

兔抗 RAB5C 多克隆抗体

中文名称：兔抗 RAB5C 多克隆抗体

英文名称：Anti-RAB5C rabbit polyclonal antibody

别名：RABL, L1880, RAB5L, RAB5CL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RAB5C

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein is a member of the Rab protein family. Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment. The function is protein transport, probably involved in vesicular traffic.
Applications:	ELISA, WB
Name of antibody:	RAB5C

Immunogen:	Synthetic peptide of human RAB5C
Full name:	Ras-related protein Rab-5C
Synonyms :	RABL, L1880, RAB5L, RAB5CL
SwissProt:	P51148
ELISA Recommended dilution:	1000-5000
WB Predicted band size:	24 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	200-1000

