

## 兔抗 RAB9A 多克隆抗体

中文名称：兔抗 RAB9A 多克隆抗体

英文名称： Anti-RAB9A rabbit polyclonal antibody

别名： RAB9

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RAB9A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Ras-related protein Rab-9A is a protein that in humans is encoded by the RAB9A gene. Rab 9A is a 201 amino acid protein that localizes to the cytoplasmic side of the cell membrane, as well as to the membrane of the Golgi apparatus and the ER, and is involved in the transport of proteins between endosomes and the trans Golgi network.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RAB9A
<b>Immunogen:</b>	Fusion protein of human RAB9A

<b>Full name:</b>	RAB9A, member RAS oncogene family
<b>Synonyms :</b>	RAB9
<b>SwissProt:</b>	P51151
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human liver cancer and human breast cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	23 kDa
<b>WB Positive control:</b>	Hela, hepg2, RAW264.7 and NIH/3T3 cells
<b>WB Recommended dilution:</b>	1000-5000



