

## 兔抗 GKAP1 多克隆抗体

中文名称：兔抗 GKAP1 多克隆抗体

英文名称： Anti-GKAP1 rabbit polyclonal antibody

别名： G kinase anchoring protein 1; FKSG21; GKAP42

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GKAP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is highly similar to the mouse cGMP-dependent protein kinase anchoring protein 42kDa. The mouse protein has been found to localize with the Golgi and recruit cGMP-dependent protein kinase I alpha to the Golgi in mouse testes. It is thought to play a role in germ cell development. Transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GKAP1
<b>Immunogen:</b>	Fusion protein of human GKAP1

<b>Full name:</b>	G kinase anchoring protein 1
<b>Synonyms:</b>	FKSG21; GKAP42
<b>SwissProt:</b>	Q5VSY0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	90-450



