

兔抗 NARFL 多克隆抗体

中文名称: 兔抗 NARFL 多克隆抗体

英文名称: Anti-NARFL rabbit polyclonal antibody

别名: PRN; HPRN; IOP1; LET1L

储存: 冷冻 (-20℃) 避光

抗原: NARFL

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

| | |
|--------------------------|---|
| Background: | Component of the cytosolic iron-sulfur protein assembly (CIA) complex, a multiprotein complex that mediates the incorporation of iron-sulfur cluster into extramitochondrial Fe/S proteins. Seems to negatively regulate the level of HIF1A expression, although this effect could be indirect. |
| Applications: | WB |
| Name of antibody: | NARFL |
| Immunogen: | Synthesized peptide derived from internal of human NARFL. |
| Full name: | nuclear prelamin A recognition factor-like |
| Synonyms : | PRN; HPRN; IOP1; LET1L |

| | |
|---------------------------------|-------------------|
| SwissProt: | Q9H6Q4 |
| WB Predicted band size: | 53 kDa |
| WB Positive control: | A549 cells lysate |
| WB Recommended dilution: | 500-3000 |

