

兔抗 RAB20 多克隆抗体

中文名称：兔抗 RAB20 多克隆抗体

英文名称： Anti-RAB20 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： RAB20

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosomes that engulf pathogens, such as S.aureus and M.tuberculosis. Plays a role in the fusion of phagosomes with lysosomes.
Applications:	ELISA, WB, IHC
Name of antibody:	RAB20
Immunogen:	Fusion protein of human RAB20
Full name:	RAB20, member RAS oncogene family
SwissProt:	Q9NX57
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-100
WB Predicted band size:	26 kDa
WB Positive control:	231 cell lysate
WB Recommended dilution:	500-2000



