

兔抗 RAB22A 多克隆抗体

中文名称：兔抗 RAB22A 多克隆抗体

英文名称：Anti-RAB22A rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：RAB22A

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|--------------------------|---|
| Background: | The protein encoded by this gene is a member of the RAB family of small GTPases. The GTP-bound form of the encoded protein has been shown to interact with early-endosomal antigen 1, and may be involved in the trafficking of and interaction between endosomal compartments. |
| Applications: | ELISA, WB, IHC |
| Name of antibody: | RAB22A |
| Immunogen: | Synthetic peptide of human RAB22A |
| Full name: | Ras-related protein Rab-22A |
| SwissProt: | Q9UL26 |

| | |
|------------------------------------|--|
| ELISA Recommended dilution: | 1000-5000 |
| IHC positive control: | Human esophagus cancer and Human lung cancer |
| IHC Recommend dilution: | 50-200 |
| WB Predicted band size: | 22 kDa |
| WB Positive control: | Mouse heart tissue |
| WB Recommended dilution: | 200-1000 |



