

兔抗 RAB25 多克隆抗体

中文名称：兔抗 RAB25 多克隆抗体

英文名称： Anti-RAB25 rabbit polyclonal antibody

别名： CATX-8; RAB11C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RAB25

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	RAB proteins, such as RAB25, are members of the RAS superfamily of small GTPases that are involved in membrane trafficking. Members of the RAB11 subfamily, including RAB25, control the return of internalized membrane-associated moieties to the cell surface.
Applications:	ELISA, IHC
Name of antibody:	RAB25
Immunogen:	Fusion protein of human RAB25

Full name:	RAB25, member RAS oncogene family
Synonyms :	CATX-8; RAB11C
SwissProt:	P57735
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human thyroid cancer
IHC Recommend dilution:	40-150

