

## 兔抗 RAB26 多克隆抗体

中文名称：兔抗 RAB26 多克隆抗体

英文名称：Anti-RAB26 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：RAB26

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a small GTP-binding protein that belongs to the largest family within the Ras superfamily. RAB26 are important regulators of vesicular fusion and trafficking. The RAB family of small G proteins regulates intercellular vesicle trafficking, including exocytosis, endocytosis, and recycling.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RAB26
<b>Immunogen:</b>	Synthetic peptide of human RAB26
<b>Full name:</b>	Ras-related protein Rab-26
<b>SwissProt:</b>	Q9ULW5

<b>IHC positive control:</b>	Human gastric cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	28 kDa
<b>WB Positive control:</b>	MCF7 and A375 cells
<b>WB Recommended dilution:</b>	200-1000



