

## 兔抗 GLCCI1 多克隆抗体

中文名称: 兔抗 GLCCI1 多克隆抗体

英文名称: Anti-GLCCI1 rabbit polyclonal antibody

别名: glucocorticoid induced 1; GCTR; GIG18; TSSN1; FAM117C

抗原: GLCCI1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein of unknown function. Expression of this gene is induced by glucocorticoids and may be an early marker for glucocorticoid-induced apoptosis. Single nucleotide polymorphisms in this gene are associated with a decreased response to inhaled glucocorticoids in asthmatic patients.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	GLCCI1
<b>Immunogen:</b>	Synthetic peptide of human GLCCI1
<b>Full name:</b>	glucocorticoid induced 1
<b>Synonyms:</b>	GCTR; GIG18; TSSN1; FAM117C

<b>SwissProt:</b>	Q86VQ1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	58 kDa
<b>WB Positive control:</b>	Human cerebella tissue lysate
<b>WB Recommended dilution:</b>	500-2000

