

## 兔抗 RAB27A 多克隆抗体

中文名称：兔抗 RAB27A 多克隆抗体

英文名称： Anti-RAB27A rabbit polyclonal antibody

别名： GS2; RAM; RAB27; HsT18676

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RAB27A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	RAB27A
<b>Immunogen:</b>	Fusion protein of human RAB27A
<b>Full name:</b>	RAB27A, member RAS oncogene family
<b>Synonyms :</b>	GS2; RAM; RAB27; HsT18676
<b>SwissProt:</b>	P51159
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	K562 cells
<b>WB Recommended dilution:</b>	500-2000



