

兔抗 NBEA 多克隆抗体

中文名称：兔抗 NBEA 多克隆抗体

英文名称： Anti-NBEA rabbit polyclonal antibody

别名： BCL8B, LYST2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NBEA

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of a large, diverse group of A-kinase anchor proteins that target the activity of protein kinase A to specific subcellular sites by binding to its type II regulatory subunits. Brain-specific expression and coat protein-like membrane recruitment of a highly similar protein in mouse suggest an involvement in neuronal post-Golgi membrane traffic. Mutations in this gene may be associated with a form of autism. This gene and its expression are frequently disrupted in patients with multiple myeloma. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Additional transcript variants m

	ay exist, but their full-length nature has not been determined.
Applications:	ELISA, IHC
Name of antibody:	NBEA
Immunogen:	Synthetic peptide of human NBEA
Full name:	neurobeachin
Synonyms :	BCL8B, LYST2
SwissProt:	Q8NFP9
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human ovarian cancer and Human cervical cancer
IHC Recommend dilution:	25-100

