

## 兔抗 TRNAU1AP 多克隆抗体

中文名称：兔抗 TRNAU1AP 多克隆抗体

英文名称：Anti-TRNAU1AP rabbit polyclonal antibody

抗原：TRNAU1AP

别名：SECP43; PRO1902; TRSPAP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclona

### 技术规格

<b>Background:</b>	Involved in the early steps of selenocysteine biosynthesis and tRNA(Sec) charging to the later steps resulting in the cotranslational incorporation of selenocysteine into selenoproteins. Stabilizes the SECISBP2, EEFSEC and tRNA(Sec) complex. May be involved in the methylation of tRNA(Sec). Enhances efficiency of selenoproteins synthesis (By similarity).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TRNAU1AP
<b>Immunogen:</b>	Fusion protein of human TRNAU1AP
<b>Full name:</b>	tRNA selenocysteine 1 associated protein 1

<b>Synonyms :</b>	SECP43; PRO1902; TRSPAP1
<b>SwissProt:</b>	Q9NX07
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human prostate cancer and Human tonsil
<b>IHC Recommend dilution:</b>	25-100

