

## 兔抗 NBN 多克隆抗体

中文名称：兔抗 NBN 多克隆抗体

英文名称： Anti-NBN rabbit polyclonal antibody

别名： ATV, NBS, P95, NBS1, AT-V1, AT-V2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NBN

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Mutations in this gene are associated with Nijmegen breakage syndrome, an autosomal recessive chromosomal instability syndrome characterized by microcephaly, growth retardation, immunodeficiency, and cancer predisposition. The encoded protein is a member of the MRE11/RAD50 double-strand break repair complex which consists of 5 proteins. This gene product is thought to be involved in DNA double-strand break repair and DNA damage-induced checkpoint activation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NBN

<b>Immunogen:</b>	Synthetic peptide of human NBN
<b>Full name:</b>	nibrin
<b>Synonyms :</b>	ATV, NBS, P95, NBS1, AT-V1, AT-V2
<b>SwissProt:</b>	O60934
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

