

兔抗 RAB39A 多克隆抗体

中文名称：兔抗 RAB39A 多克隆抗体

英文名称： Anti-RAB39A rabbit polyclonal antibody

别名： RAB39A, member RAS oncogene family; RAB39

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RAB39A

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Plays a role in the maturation and acidification of phagosomes that engulf pathogens, such as S.aureus and M.tuberculosis. Plays a role in vesicular trafficking. Plays a role in the fusion of phagosomes with lysosomes. Negatively regulates LPS-induced autophagosome formation in macrophages possibly by implicating PI3K. May be involved in multiple neurite formation.
Applications:	ELISA, WB, IHC
Name of antibody:	RAB39A
Immunogen:	Fusion protein of human RAB39A

Full name:	RAB39A, member RAS oncogene family
Synonyms:	RAB39
SwissProt:	Q14964
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	25 kDa
WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	200-1000



