

## 兔抗 TRPA1 多克隆抗体

中文名称: 兔抗 TRPA1 多克隆抗体

英文名称: Anti-TRPA1 rabbit polyclonal antibody

别名: FEPS; ANKTM1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TRPA1

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The structure of the protein encoded by this gene is highly related to both the protein ankyrin and transmembrane proteins. The specific function of this protein has not yet been determined; however, studies indicate the function may involve a role in signal transduction and growth control.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TRPA1
<b>Immunogen:</b>	Synthetic peptide of human TRPA1
<b>Full name:</b>	transient receptor potential cation channel, subfamily A, member 1
<b>Synonyms:</b>	FEPS; ANKTM1

<b>SwissProt:</b>	O75762
<b>IHC positive control:</b>	Human gastric cancer and human lung cancer
<b>IHC Recommend dilution:</b>	25-100

