

兔抗 RAB40B 多克隆抗体

中文名称：兔抗 RAB40B 多克隆抗体

英文名称： Anti-RAB40B rabbit polyclonal antibody

别名： RAR; SEC4L

抗原： RAB40B

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene has similarity to a yeast protein which suggests a role of the gene product in regulating secretory vesicles.
Applications:	WB
Name of antibody:	RAB40B
Immunogen:	Synthesized peptide derived from internal of human RAB40B.
Full name:	RAB40B, member RAS oncogene family
Synonyms :	RAR; SEC4L
SwissProt:	Q12829
WB Predicted band size:	31 kDa

WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

