

兔抗 RAB41 多克隆抗体

中文名称：兔抗 RAB41 多克隆抗体

英文名称： Anti-RAB41 rabbit polyclonal antibody

别名： RAB41, member RAS oncogene family

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RAB41

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a small GTP-binding protein that belongs to the largest family within the Ras superfamily. These proteins function as regulators of membrane trafficking. They cycle between inactive GDP-bound and activated GTP-bound states, which is controlled by GTP hydrolysis-activating proteins (GAPs). This family member can be activated by the GAP protein RN-Tre, and it is localized to the Golgi complex.
Applications:	ELISA, WB, IHC
Name of antibody:	RAB41

Immunogen:	Fusion protein of human RAB41
Full name:	RAB41, member RAS oncogene family
SwissProt:	Q5JT25
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human thyroid cancer
IHC Recommend dilution:	100-200
WB Predicted band size:	25 kDa
WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	500-2000



