

## 兔抗 RABEP2 多克隆抗体

中文名称: 兔抗 RABEP2 多克隆抗体

英文名称: Anti-RABEP2 rabbit polyclonal antibody

别名: FRA

抗原: RABEP2

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human Mouse Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	Plays a role in membrane trafficking and in homotypic early endosome fusion.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	RABEP2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human RABEP2.
<b>Full name:</b>	rabaptin, RAB GTPase binding effector protein 2
<b>Synonyms :</b>	FRA
<b>SwissProt:</b>	Q9H5N1

<b>WB Predicted band size:</b>	64 kDa
<b>WB Positive control:</b>	Hela cells, Jurkat cells, 293 cells and H uvec cells lysates
<b>WB Recommended dilution:</b>	500-3000

