

## 兔抗 NCAPG2 多克隆抗体

中文名称：兔抗 NCAPG2 多克隆抗体

英文名称：Anti-NCAPG2 rabbit polyclonal antibody

别名：MTB; CAPG2; LUZP5; CAP-G2; hCAP-G2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NCAPG2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that belongs to the Condensin2nSMC family of proteins. The encoded protein is a regulatory subunit of the condensin II complex which, along with the condensin I complex, plays a role in chromosome assembly and segregation during mitosis. A similar protein in mouse is required for early development of the embryo. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NCAPG2
<b>Immunogen:</b>	Synthetic peptide of human NCAPG2

<b>Full name:</b>	non-SMC condensin II complex, subunit G2
<b>Synonyms :</b>	MTB; CAPG2; LUZP5; CAP-G2; hCAP-G2
<b>SwissProt:</b>	Q86XI2
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	50-200

