

兔抗 NCAPG 多克隆抗体

中文名称：兔抗 NCAPG 多克隆抗体

英文名称： Anti-NCAPG rabbit polyclonal antibody

别名： non-SMC condensin I complex subunit G; CAPG; CHCG; YCG1; NY-MEL-3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NCAPG

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a subunit of the condensin complex, which is responsible for the condensation and stabilization of chromosomes during mitosis and meiosis. Phosphorylation of the encoded protein activates the condensin complex. There are pseudogenes for this gene on chromosomes 8 and 15. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	NCAPG
Immunogen:	Fusion protein of human NCAPG
Full name:	non-SMC condensin I complex subunit G

Synonyms:	CAPG; CHCG; YCG1; NY-MEL-3
SwissProt:	Q9BPX3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human ovarian cancer
IHC Recommend dilution:	30-150



