

## 兔抗 TRPM8 多克隆抗体

中文名称：兔抗 TRPM8 多克隆抗体

英文名称：Anti-TRPM8 rabbit polyclonal antibody

别名：TRPP8; LTRPC6

相关类别：一抗

储存：冷冻 (-20℃)

抗原：TRPM8

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Transient receptor potential ion channels (TRPC) are a superfamily of six transmembrane segment-spanning gated cation channels. TRP subtypes mediate store-operated Ca <sup>2+</sup> entry, a process involving Ca <sup>2+</sup> influx and replenishment of Ca <sup>2+</sup> stores formerly emptied through the action of inositol 1,4,5-trisphosphate production and other Ca <sup>2+</sup> mobilizing agents.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TRPM8
<b>Immunogen:</b>	Fusion protein of human TRPM8
<b>Full name:</b>	transient receptor potential cation channel, subfamily

	y M, member 8
<b>Synonyms :</b>	TRPP8; LTRPC6
<b>SwissProt:</b>	Q7Z2W7
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

