

兔抗 NCAPH 多克隆抗体

中文名称: 兔抗 NCAPH 多克隆抗体

英文名称: Anti-NCAPH rabbit polyclonal antibody

别名: BRRN1; CAP-H

储存: 冷冻 (-20℃)

抗原: NCAPH

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the barr gene family and a regulatory subunit of the condensin complex. This complex is required for the conversion of interphase chromatin into condensed chromosomes. The protein encoded by this gene is associated with mitotic chromosomes, except during the early phase of chromosome condensation. During interphase, the protein has a distinct punctate nuclear localization. Alternatively spliced transcript variants encoding different proteins have been described.
Applications:	ELISA, IHC
Name of antibody:	NCAPH

Immunogen:	Fusion protein of human NCAPH
Full name:	non-SMC condensin I complex subunit H
Synonyms:	BRRN1; CAP-H
SwissProt:	Q15003
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100



