

兔抗 Biotin 多克隆抗体

- 中文名称: 兔抗 Biotin 多克隆抗体
- 英文名称: Anti-Biotin rabbit polyclonal antibody
- 别 名: Biotin; Vitamin H; Vitamin B7
- 相关类别: 标签抗体
- 储 存: 冷冻 (-20℃)
- 宿 主: Rabbit
- 抗 原: Biotin
- 反应种属: Biotin
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the i nitial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is req uired by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes invol ved in incorporation of carbon dioxide into organic co mpounds. Biotin has a MW of 244 Da.
Applications:	ELISA
Name of antibody:	Biotin



Immunogen:	KLH conjugated Biotin
Full name:	Biotin
ELISA Recommended dilution:	5000-10000